

Corning® 88-581-CM Medium for Activation and Expansion of Human CIK Cells

Protocol

CORNING

Corning 88-581-CM medium is part of the Corning 500 media series, which is widely used in adoptive immunotherapy. It is specifically designed for activation and expansion of cytokine induced killer (CIK) cells. 88-581-CM is manufactured with high quality reagents and GMP-grade raw materials. No other proteins are present in the medium except the injectable level of serum albumin and recombinant human insulin.

Content and Storage

Description	Size	Qty/Pk	Shelf Life
88-581-CM	1,000 mL	1	12 months

Reagents and Materials

- ▶ Lymphocyte Separation Medium (LSM, Corning Cat. No. 25-072-CI)
- ▶ PBS (Corning Cat. No. 21-040-CV)
- ▶ INF- γ (Peprotech Cat. No. 300-02)
- ▶ IL-2 (Corning Cat. No. 354043)
- ▶ Anti-CD3 antibody (OKT3 provided by KOHJIN BIO)
- ▶ T-75 flask, TC-treated (Corning Cat. No. 430641)
- ▶ T-225 flask, TC-treated (Corning Cat. No. 431081)
- ▶ Gas-permeable cell culture bag (Corning Cat. No. 88-610-20)

CIK Cells Activation and Expansion

1. On Day 0, centrifuge the anti-coagulated blood at 400 x g for 10 minutes at room temperature (RT) and then transfer the plasma (on the top layer) into a new tube. Inactivate the auto-plasma at 56°C for 30 minutes and then centrifuge at 800 x g for 20 minutes. Store the supernatant at 4°C for further use. Add equal volume of PBS (without calcium and magnesium) as a replacement of the removed auto-plasma to maintain a constant volume and resuspend the haemocytes gently.
Note: Pre-adding 0.1% human serum albumin (HSA) in PBS is helpful to maintain haemocytes viability (optional step).
2. Prepare peripheral blood mononuclear cells (PBMCs) from the above blood sample using LSM according to the manufacturer's directions.
Note: Use freshly collected human blood (within 2 hours of collection) for better performance; do not use blood that is drawn more than 24 hours prior to use.
3. Wash the PBMCs once with at least 5-fold PBS buffer and then centrifuge the sample at 500 x g for 10 minutes at RT.
Note: Pre-adding 0.1% HSA in PBS is helpful to maintain PBMCs viability (optional step).
4. Adjust the density of PBMCs to 1 to 2 x 10⁶ cells/mL using 88-581-CM medium containing 1,000 IU/mL INF- γ and 5% auto-plasma (from Step 1).
Note: Normally, more than 1 x 10⁶ lymphocytes can be obtained from 1 mL blood sample.
5. Seed cells into desired culture vessel and then incubate the vessel at 37°C in a humidified atmosphere of 5% CO₂ in air.

6. After 24 hours, add IL-2 and anti-CD3 antibody at final concentrations of 500 IU/mL and 50ng/mL, respectively.
7. Around Day 5, adjust the cell density to 1×10^6 cells/mL by using 88-581-CM medium containing 500 IU/mL IL-2 and 5% auto-plasma.
8. Based on the cell proliferation status (typically on Day 7 before cells enter the exponential growth phase), transfer the cell suspension to a larger flask or gas-permeable culture bag for large-scale cell expansion. Use fresh 88-581-CM medium containing 500 IU/mL IL-2 and 0.5% to 1% auto-plasma to dilute the cell suspension.
Note: Ensure that the final cell density is more than 5×10^5 cells/mL.
9. Depending on the cell proliferation status, replenish fresh medium containing 500 IU/mL IL-2 every 2 to 3 days; addition of 1% auto-plasma is recommended until all the auto-plasma has been exhausted.
Note: Keep the final cell density within 1 to 2×10^6 cells/mL after each dilution.
10. Harvest mature CIK cells around Day 14.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit www.corning.com/lifesciences or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

CORNING

Corning Incorporated
Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

ASIA/PACIFIC

Australia/New Zealand
t 61 427286832

China
t 86 21 3338 4338
f 86 21 3338 4300

India
t 91 124 4604000
f 91 124 4604099

Japan

t 81 3-3586 1996
f 81 3-3586 1291

Korea

t 82 2-796-9500
f 82 2-796-9300

Singapore

t 65 6572-9740
f 65 6861-2913

Taiwan

t 886 2-2716-0338
f 886 2-2516-7500

EUROPE

CSEurope@corning.com

France

t 0800 916 882
f 0800 918 636

Germany

t 0800 101 1153
f 0800 101 2427

The Netherlands

t 020 655 79 28
f 020 659 76 73

United Kingdom

t 0800 376 8660
f 0800 279 1117

All Other European Countries

t +31 (0) 206 59 60 51
f +31 (0) 206 59 76 73

LATIN AMERICA

grupoLA@corning.com

Brasil

t 55 (11) 3089-7400

Mexico

t (52-81) 8158-8400