User Manual

(D) DIGI-SENSE. Digital Pocket Thermometer

with NIST-Traceable Calibration

Model 20250-32



Introduction

The Digi-Sense Large-Head Pocket Thermometer (Model 20250-32) offers fast response and laboratory accuracy with 0.1°F (0.1°C) resolution over the entire measurement range. Advanced features include °C/°F unit selection and automatic power-off. The instrument is fully tested and calibrated to NIST-traceable standards. Careful use of this meter will provide years of reliable service.

Unpacking

Check individual parts against the list of items below. If anything is missing or damaged, please contact your instrument supplier immediately.

- 1. Pocket thermometer
- 2. One LR44 (1.5V) battery, installed
- 3. Protective sleeve
- 4. User manual
- 5. NIST-traceable calibration report with data

Key Features

- User-selectable °C or °F units
- · Automatic power-off to conserve battery life

Meter Description

- 1. Protective cover
- 2. Temperature probe
- 3. LCD
- 4. °C/°F button
- 5. Power ON/OFF button



Setup and Operation

Power

- 1. Remove the thermometer from the protective sleeve.
- 2. Press the Power on/off button to activate.

Taking Measurements

Insert the stem of the thermometer at least $\frac{1}{2}$ " into the testing area.

Note: When the thermometer is initially powered on, the LCD will display "- 0.0 °C" or "32.0 °F" for two seconds. Then the actual temperature will display. To save battery life, the thermometer will automatically turn off after approximately 20 minutes.

Temperature Measurements

Press the °C/°F button to change the temperature unit from °C to °F.

Specifications

Temperature range	Full usable range: -40 to 482°F (-40 to 250°C) NIST certified range: 32 to 482°F (0 to 150°C)
Resolution	0.1°F (0.1°C)
Accuracy	1.5% ±4°F (±2°C) full scale

Display LCD

Response time 20 seconds

Unit shuts off automatically after 20 minutes Auto power-off

to preserve battery life

Weight 1.5 oz (42.5 a)

Dimensions (L x W x D) 6" x 11/2" x 11/4" (150 x 40 x 28 mm)

One LR44 (1.5V) battery Power

Maintenance, Recalibration, and Repair

Cleaning and Storage

- The thermometer should be cleaned with a damp cloth and mild detergent when necessary. Do not use solvents or abrasives
- Store the meter in an area with moderate temperature and humidity.

Battery Replacement

When the battery power falls low, open the battery compartment on the rear of the unit and replace the LR44 battery. Close the cover securely when finished. It is recommended that Digi-Sense products are calibrated annually to ensure proper function and accurate measurements; however, your quality system or regulatory body may require more frequent calibrations. To schedule your recalibration, please contact InnoCal, an ISO 17025 calibration laboratory accredited by A2LA.



Phone: 1-866-INNOCAL (1-866-466-6225)

Fax: 1-847-327-2993

E-mail: sales@innocalsolutions.com

Web: InnoCalSolutions.com

For Product and Ordering Information, Contact:



Toll-Free: 1-800-323-4340 Phone: 1-847-549-7600 Fax: 1-847-247-2929

ColeParmer.com/Digi-Sense



Toll-Free: 1-800-358-5525 Phone: 1-847-327-2000 Fax: 1-847-327-2700

Davis.com/Digi-Sense