



Cultureware

Sterile, Single-Use Erlenmeyer Flasks

Choose the size that fits your needs

The Thermo Scientific™ Sterile, Single-Use Erlenmeyer Flasks with Plain Bottom are ideal for suspension of cell culture, and media preparation, mixing, and storage. Made of crystal-clear polyethylene terephthalate glycol (PETG), these flasks are robust and break resistant. The vented cap is designed for continuous gas exchange while maintaining sterility and preventing leakage.

Product features include:

- Molded-in graduation for accurate volume assessment
- Vented cap to allow gas exchange
- Sterilized by gamma irradiation per ISO 11137-2
- Non-pyrogenic and non-cytotoxic
- DNase-, RNase-, and human DNA-free
- BSE and TSE risk-free

Ordering information

Capacity	Outer diameter	Case size	Cat. No.
125 mL	68 mm	24	4117-0125
250 mL	90 mm	12	4117-0250
500 mL	100 mm	12	4117-0500

Buy now

Learn more about the entire culture flask portfolio at thermofisher.com/flasks

thermo scientific