










SECTION C-C



SECTION A-A




SECTION D-D



DETAIL B
SCALE 10:1

ISSUE	DATE	AMENDMENT
A	25.04.2016	Original drawing


www.4ti.co.uk


DO NOT SCALE
IF IN DOUBT ASK

CUSTOMER:	
CUSTOMER DRG. No.:	ISSUE:
MATERIAL: FRAME: Polycarbonate TUBES: Polypropylene	
TOLERANCES:	0 ± 0.25 0.0 ± 0.20 0.00 ± 0.10 ANGULAR ± 0°30'
GENERAL TOLERANCES (UNLESS OTHERWISE STATED)	(ALL DIMENSIONS IN MILLIMETRES)

DRAWN:	ZE
DATE:	25.04.2016
CHECKED:	HO
DATE:	
SCALE:	

TITLE: FrameStar® Break2Way Low Profile Two direction grooves Marketing drawing	
DRAWING No.:	4ti-1400
SHEET No.:	1 OF 1
ISSUE:	A