





SECTION A-A



TOOL INSERT TO
ALLOW CUSTOMERS
NAME TO BE ENGRAVED
IF REQUESTED





SECTION B-B



DETAIL C
SCALE 10 : 1

ISSUE	DATE	AMENDMENT

 www.4ti.co.uk	 DO NOT SCALE IF IN DOUBT ASK	CUSTOMER:	SIGNATURE:	DRAWN:	KIRAN. S	TITLE: 0291-LID Marketing Drawing		
		CUSTOMER DRG. No.:	ISSUE:	DATE:	24/05/2019			
		MATERIAL:	GLASS FILLED POLYCARBONATE		CHECKED:	HANNA. O		
		TOLERANCES:			DATE:		DRAWING No.:	
		GENERAL TOLERANCES (UNLESS OTHERWISE STATED)			SCALE:		SHEET No.:	
0 ±0.15 0.0 ±0.05 0.00 ±0.02 ANGULAR ±0.30° (ALL DIMENSIONS IN MILLIMETRES)					PO-0291-LID			
					1 OF 1			
					ISSUE: B			