

Clifton Stirred Water Baths

NE4 - D/CT Series Stirred Digital Baths

Clifton Range®



NE4-22D/CT



Clifton NE4D/CT digital stirred baths come complete with sturdy clear polycarbonate tank, allowing observation of samples. Complete with stainless steel shelf designed for good water circulation and temperature uniformity. All Clifton stirred baths feature a low liquid level float switch, which ensures heating is switched off if the liquid falls below the recommended level.

NE4 - 14D/CT

NE4 - 22D/CT

GENERAL			
Temp. Range	Ambient + 5°C - 90°C (0°C with accessory cooling)		
Temp. Alarm	Under/over temperature alarm is automatically set 4°C above and 4°C below set point		
Temp. Control	PID Digital, LED display, resolution 0.1°C		
Sensitivity/Uniformity	±0.05°C / ±0.01°C		
Timer	Visual indicator and audible end message - Variable 00:01 to 99:59 hh:mm		
Heater Power (220-240V)	1250W		
Voltage	120V or 230V models available 50-60Hz		
DIMENSIONS			
Tank Capacity	14 Litre	22 Litre	
External (mm)	530w x 325d x 295h	530w x 325d x 345h	
Internal (mm)	385w x 290d x 150h	385w x 290d x 200h	
Water Level Max	130		

NE4-D/CT Series Stirred Digital Baths

Accessories



Polypropylene Spheres (Pack of 200)

Bath	NE4-14D/CT	NE4-22D/CT
No. of Packs	2 Packs	3 Packs
Accessory Code	BP0368	



Immersion Dip Cooler [ambient to -20°C] *

Bath	NE4-14D/CT	NE4-22D/CT
Accessory Code	DC1-300	

* Please note : NE4D/CT range will only operate to 0°C with accessory cooling.



Stainless Steel Test Tube Racks

Bath	NE4-14D/CT	NE4-22D/CT
No. of Racks	3	5
Accessory Code	6870 : 26 holes x 17mm DIA	
	6871 : 16 holes x 26mm DIA	
	6872 : 36 holes x 13mm DIA	
	6873 : 18 holes x 19mm DIA	
	6900 : 12 holes x 32mm DIA	

Universal Accessories

CAL	Calibration Record
BX0616	Draining Syphon