

# Clifton Stirred Water Baths NE4 - T Series Analogue Thermostatic Stirred Baths

*Clifton Range*<sup>®</sup>



NE4-14T

Clifton NE4-T thermostatically controlled water baths are designed for immersion of flasks/ bottles/racks in a stable temperature environment. The stainless steel tank is contained within a sturdy outer case and supplied complete with a stainless steel shelf, which provides good water circulation and temperature stability. All Clifton stirred baths feature a low liquid level float switch, which ensures heating is switched off if the liquid falls below the recommended level.

## NE4 - 8 T

## NE4 - 14T

GENERAL	
Temp. Range	Ambient + 5°C - 99°C
Temp. Control	Analogue dial pointer, graduated in 5°C divisions
Sensitivity/Uniformity	±0.2°C / ±0.05°C
Heater Power (220-240V)	1250W
Voltage	230V models available 50-60Hz
DIMENSIONS	
Tank Capacity	8 Litre
External (mm)	271w x 332d x 170h
Internal (mm)	129w x 300d x 150h
Water Level Max. (mm)	130

## NE4 - 22T

## NE4 - 28T

GENERAL	
Temp. Range	Ambient + 5°C - 99°C
Temp. Control	Analogue dial pointer, graduated in 5°C divisions
Sensitivity/Uniformity	±0.2°C / ±0.05°C
Heater Power (220-240V)	1250W
Voltage	230V models available 50-60Hz
Drain	N/A
DIMENSIONS	
Tank Capacity	22 Litre
External (mm)	537w x 332d x 170h
Internal (mm)	395w x 298d x 150h
Water Level Max. (mm)	130

# NE4-T Series Analogue Thermostatic Stirred Baths Accessories



## Stainless Steel Gabled lid with Insulated Handle

Bath	NE4-8T	NE4-14T	NE4-22T	NE4-28T
Accessory Code	SL4-8	SL4-14	SL4-22	



## Stainless Steel Flat One Piece lid with Insulated Handle

Bath	NE4-8T	NE4-14T	NE4-22T	NE4-28T
Accessory Code	LD4-8	LD4-14	LD4-22	



## Polypropylene Spheres (Pack of 200)

Bath	NE4-8T	NE4-14T	NE4-22T	NE4-28T
No. of packs	1 Pack	2 Packs	3 Packs	
Accessory Code	BP0368			



## Stainless Steel Raised Shelves

Bath	NE4-8T	NE4-14T	NE4-22T	NE4-28T
Accessory Code	N/A	RS4-14	RS4-22	



## Drain Tap

Bath	NE4-8T	NE4-14T	NE4-22T	NE4-28T
Accessory Code	N/A			BX0953

# NE4-T Series Analogue Thermostatic Stirred Baths

## Accessories



### Stainless Steel Test Tube Racks

Bath	NE4-8T	NE4-14T	NE4-22T	NE4-28T
No. of Racks	1	3	5	
Accessory Code	6870 : 26 holes x 17mm DIA			
	6871 : 16 holes x 26mm DIA			
	6872 : 36 holes x 13mm DIA			
	6873 : 18 holes x 19mm DIA			
	6900 : 12 holes x 32mm DIA			

## Universal Accessories

BX0616	Draining Syphon
--------	-----------------