

THE NEW **DURAN®** GL 45 BOTTLE CARRYING SYSTEM



DURAN[®] GL 45 BOTTLE CARRYING **SYSTEM**

The new modular Bottle Carrying System fits standard laboratory bottles that have a GL 45 neck finish and is compatible with four bottle sizes 2, 5, 10 and 20 litres. It is ergonomically designed, and offers health & safety benefits to the user and allows easy and safe transport of bottles around the laboratory or production facility. It makes pouring liquids from the larger bottles more comfortable and safer. All materials are fully autoclavable.



GL 45 Bottle Collar

- The collar is securely attached by pushing down over the GL 45 thread with 4 teeth that lock under transfer ring.
- The collar is only removable by cutting so attachment is 100% secure in use. A cut line indicator for removal of collar when required.

Horizontal Straps

- Durable stitched construction
- Made from "anti-slip" material



Silicone Support Base

- larger 5, 10 and 20 litre bottles
- Protection of glass bottle base from
- Spacer to protect glass from
- Moulded-in holes minimize liquid

Lot Specific Retrace Code

for traceability and certification



Optional Extra: Metal Dolly

- Stainless steel dolly for 10 or 20 litre bottles (2 sizes only)
- Fully compatible with carrying system
- Easy to clean, crevice-free construction

FOR EASY AND SAFE **CARRYING & POURING**



Easy and safe carrying using handle

- Handle attaches at the neck and fully
- supports the base
- For 2, 5, 10 and 20 litre bottles



Safe and controlled pouring using handle • Safer pouring of liquids using a second lower "Pouring Handle" located on the side strap

• Available on the 5, 10 and 20 litre carriers only

Ordering information

| Description | Pack Unit | Cat. No. |
|---|-----------|-----------|
| DURAN® GL 45 Bottle Carrying System for the 2 Litre DURAN® Bottle (Blue) | 1 | 292836307 |
| DURAN® GL 45 Bottle Carrying System for the 2 Litre DURAN® Bottle (Purple) | 1 | 292826306 |
| DURAN® GL 45 Bottle Carrying System for the 2 Litre DURAN® Bottle (Red) | 1 | 292816305 |
| DURAN® GL 45 Bottle Carrying System for the 2 Litre DURAN® Bottle (Yellow) | 1 | 292806304 |
| DURAN® GL 45 Bottle Carrying System for the 5 Litre DURAN® Bottle (Yellow) | 1 | 292807309 |
| DURAN® GL 45 Bottle Carrying System for the 10 Litre DURAN® Bottle (Yellow) | 1 | 292808605 |
| DURAN® GL 45 Bottle Carrying System for the 20 Litre DURAN® Bottle (Yellow) | 1 | 292809104 |
| DURAN® Metal Dolly for 10 L DURAN® Bottle (stainless steel) available July 2018 | 1 | 299018607 |
| DURAN® Metal Dolly for 20 L DURAN® Bottle (stainless steel) available July 2018 | 1 | 299019106 |

DURAN[®] Metal Dolly for 20 L DURAN[®] Bottle (stainless steel) available July 2018

* Bottles supplied separately

Fast and secure tightening

of the polypropylene straps with two side strap adjusters for the larger 5, 10 and 20 litre bottles

> 4 sizes available: 2 L / 5 L / 10 L / 20 L



2 L bottle carrier in 4 colors for color coding









DWK Life Sciences GmbH Hattenbergstraße 10 55122 Mainz, Germany

3

3

Phone: +49 6131 - 1445 4131 Fax: +49 6131 - 1445 4016 sales@DWK-LifeSciences.com www.DWK-LifeSciences.com



ETRACE