


PASTEUR PIPETTES
MADE OF
BOROSILICATE GLASS 5.1




DURAN
WHEATON
KIMBLE

Excellence in your hands

PASTEUR PIPETTES MADE OF BOROSILICATE GLASS 5.1

Due their material properties, Pasteur pipettes made of borosilicate glass 5.1 are a high-quality alternative to plastic pipettes and pipettes made of soda-lime glass.



Applications

- Pasteur pipettes are intended for one-time use in microbiological and serological laboratories
- Ideal for distributing aliquots quantities of liquid

Features & Benefits

- Glass pasteur pipettes with short and long drawn fine tip
- Made of clear borosilicate glass 5.1
- USP Type I neutral glass
- Pipettes conform to ISO 7712
- Inured to rapid temperature changes and chemical attack
- Very good chemical resistance to acidic and alkaline solutions
- Borosilicate glass 5.1 is more stable than soda-lime glass and therefore ideal for the transfer of hot liquids
- Virtually inert borosilicate glass minimizes interaction with the material and reduces the risk of sample contamination
- Rubber teat, made from natural rubber, can be ordered separately

Ordering information

	Description	l (mm)	Pipette d (mm)	Tip d (mm)	Volume (mL)	Pack Unit	Cat. No.
NEW	Pasteur pipette	150	6.95	approx. 1.1	approx. 2	1000 = 4 x 250	23 330 15 04
NEW	Pasteur pipette	230	6.95	approx. 1.1	approx. 2	1000 = 4 x 250	23 330 23 03

Accessory parts

Description	d (mm)	h (mm)	Pack Unit	Cat. No.
Rubber teat made from natural rubber	15	35	100	29 200 01 02



DWK Life Sciences GmbH
Hattenbergstrasse 10
55122 Mainz
Germany

Phone: +49 6131 - 1445 4131
Fax: +49 6131 - 1445 4016
sales@DWK-LifeSciences.com
www.DWK-LifeSciences.com

