

PRODUCT INFORMATION

Dual Approved Powderfree Nitrile Gloves

Description:

When strength and dexterity are of the utmost importance, our allergy-free Powderfree Nitrile Gloves are the ideal choice. From food handling to engineering, these versatile gloves are sure to suit the task at hand.

- ✓ Medical Grade AQL 1.5
- ✓ Conforms to EN 455
- ✓ BS EN 374-1/Type C
- ✓ Tested in accordance with BS EN 420:2003+A1:2009
- ✓ Latex free
- ✓ Dual approved for medical and industrial use
- ✓ Excellent barrier protection
- ✓ Textured finger tips
- ✓ Moulds to shape of hand
- ✓ Can be worn under working gloves
- ✓ Low chemical contact use
- ✓ Protection against bacteria and fungi
- ✓ Protection against viruses

Type:

Nitrile
Powder free

Applications:

Clinical tasks, laboratories, food processing, food handling, automotive and engineering

Packaging: 100 pieces per dispenser - 10 dispensers per case



	EN ISO 374-1:2016/Type C	Chemical Class	Levels	EN 374-4:2013 Degradation
K	40% Sodium Hydroxide	Inorganic Base	6	3.5%

Product Codes:

12611-4 S - XL Nitrile Powder Free - Blue

Colourways:

 Blue

PPE Regulation (EU) 2016/425


Plastic Materials and Articles in Contact with Food Regulation (EU) No 10/2011

EN 374-1:2016/Type C

EN 374-4:2013 % Degradation

EN 374-5:2016

Standards:



Size Chart:

Measurements in mm	7 (S)	8 (M)	9 (L)	10 (XL)	11 (2XL)
Hand Length (A)	240	240	240	240	240
Hand Circumference (B)	178	203	229	254	279

Issue No.: 5

Issue Date: 02.12.2021