

JB Nova - 5L, 12L, 18L and 26L

NEW

PERFORMANCE

› Set and Forget™ technology - fast heat-up, reliable temperature control

INTUITIVE

> Simple to use controls, clear bright display

QUALITY

3 year warranty, UK design and manufacture

PRACTICAL

> Display lock - disables front panel controls

CONVENIENT

> Drain tap, tray and lid included. Suitable for use with heat transfer beads

VALUE

> Advanced dry start and run dry protection - prevents costly service repairs







Applications

- Education school/college university teaching lab
- Healthcare sample warming
- Industry QA/QC testing

Ideal choice for schools and colleges

Specifications

	JBN5	JBN12	JBN18	JBN26
C€	3 kg l: 215 mm w: 335 mm h: 200 mm	6 kg I: 380 mm w: 360 mm h: 225 mm	9 kg I: 590 mm w: 335 mm h: 275 mm	9 kg I: 590 mm w: 335 mm h: 275 mm
Capacity	5L	12L	18L	26L
Temperature range	ambient +5 to 95°C			
Temperature display and setting resolution	0.5°C			
Stability (DIN 12876-3 @70 °C)	±0.5°C			
Working area (I x w) mm	131 x 281	281 x 306	485 x 281	481 x 278
User defined calibration points	1	1	1	1
Drain tap	-	•	•	•
Heater power 230V/120V (kW)	0.35 / 0.35	0.8 / 0.8	1.4 / 1.05	1.4 / 1.05
Supply voltage (V)	120 or 230 (50-60Hz)			

Accessories (full range available on website)

Flat lid with ring sets

Variable hole diameter to accommodate tall vessels whilst reducing evaporation



Polypropylene spheres

Useful alternative to a lid, minimises evaporation and heat loss. Easy access to vessels in the bath; particularly useful for tall vessels



Test tube racks

Stainless steel - range of sizes available for 0.5 to 1.5 ml microtubes or 10 to 30 mm \varnothing tubes



Raised shelves

Stainless steel - shelf covers half area of bath. Reversible, allows two shelf depths



Drain stopper

Recommended when using heat transfer beads, to prevent beads entering drain hole



Alternative water bath ranges:



JB Nova - basic

- Simple and practical, great value
- 3 bath sizes: 5L, 12L, 18L

SUB Aqua Pro - advanced range

- Temperature range amb. +5 to 99°C
- Stability ±0.2°C
- 8 bath sizes: 2L, 2L shallow, 5L, 12L, 18L, 26L, 5L & 12L dual bath
- 3 temperature presets
- User defined high temperature alarm and cutout



Call us:

+ 44 (0) 1763 264 741

Email:

salesdesk@grantinstruments.com