

# General Purpose pH Electrode | HI-1217D

## Ideal for general use

The HI-1217D is a general purpose pH electrode.

This electrode has a 1m cable and DIN connector.

- diam: 12mm
- length: 110mm
- temp: 0 to 70°C
- max pressure: 2 bar

### Recommended use

- general purpose

### Ordering Information


**General Purpose pH Electrode | HI-1217D**  
supplied with DIN connector and 1m cable

#### Electrode Storage Solutions

electrode storage solution (500 mL) HI-70300L  
electrode storage solution (230 mL) HI-70300M

### Specifications

Product Code	HI-1217D
Reference	single, Ag/AgCl
Junction	ceramic, single
Connector	DIN
Electrolyte	gel
Max Pressure	2 bar
Range	pH: 0 to 13; 0 to 70°C
Tip Shape	spheric (dia: 5.0 mm)
Temperature Sensor	yes
Amplifier	yes
Body Material	PEI
Recommended Use	general purpose



**DIN  
Connection**



Authorised distributor



Scientific Laboratory Supplies Ltd.  
Orchard House  
The Square, Hessle  
East Yorkshire, HU13 0AE

Tel: +44 1482 496 666  
Fax: +44 1482 649 667  
Email: slsexport@scientific-labs.com  
Web: www.scientific-labs.com