




Economical, Reusable & Moisture Resistant Plastic

True north® freezer boxes by Heathrow Scientific, were designed with sustainability in mind and are made from a specially formulated corrugated or thin film polypropylene. True north freezer boxes are moisture resistant, which significantly reduces the chance of mold growth. With superior durability, polypropylene outlasts the cardboard box alternatives that are susceptible to softening and producing dust and can freeze (down to temperatures of -80°C) and thaw multiple times without a problem.

When it comes to saving valuable lab space, true north® Flatpack Freezer Boxes are supplied “flat-packed”, which eliminates the hassle of stacking and reduces the amount of space needed for storage. When ready to use, assembly is easy and pre-assembled dividers are included with each box. True north freezer boxes come complete in packs of 10, and multiple sizes and colors available.



Freezer boxes are supplied “flat-packed”
Eliminating the hassle of stacking and reduces the amount
of space needed for storage



