

droplet™

Lightweight, roomy, and easy-to-clean blood collection tray

Blood Collection Tray, Inserts and Lid

Keeping your tray clean and organized has never been easier. Along with the sturdy plastic frame you get a set of economical, disposable inserts specifically designed for blood collection. The replaceable tray insert comes in two configurations and includes an extra-deep compartment for the sharps container.

Other tray features include a supply drawer, side pockets, and a ratcheting handle with soft rubber grip that rotates out of the way for easy access to the tray contents.

Each tray comes with its own 13mm or 16mm Heathrow OneRack™ test tube rack (your choice) and two disposable polystyrene inserts. Additional inserts sold in packs of four.


Provide extra protection for the contents of the Droplet Tray with the Droplet lid. The custom designed sturdy lid offers a snug fit, protects contents and provides access to the tray handle in its vertical position. The lids are stackable, made from polystyrene, and cable ties can be used for extra security.

Droplet™ Configurations


Part No.	Description	L x W x H	Handle Height
HS2200A	Droplet™ blood collection tray with 2 inserts (style a) and a 36-place, 13 mm test tube rack	16.8 x 10.4 x 6.7 in (42.9 x 26.5 x 17.2 cm)	11 in (28 cm)
HS2200B	Droplet™ blood collection tray with 2 inserts (style b) and a 36-place, 16 mm test tube rack	16.8 x 10.4 x 6.7 in (42.9 x 26.5 x 17.2 cm)	11 in (28 cm)

Droplet™ Inserts

Part No.	Description	L x W
HS2201A	Droplet™ blood collection tray insert A for 13 mm tube rack	14.2 x 8.2 in (36 x 21 cm)
HS2201B	Droplet™ blood collection tray insert B for 16 mm tube rack	14.2 x 8.2 in (36 x 21 cm)



Disposable insert A
for 13 mm tube rack



Disposable insert B
for 16 mm tube rack



Handle with soft rubber grip for easy, safe transport



Droplet™ Lid

Cover for the Droplet™ Blood Collection Tray

Protect the contents of the Droplet Tray with this Droplet lid. Custom designed sturdy lid is a snug fit, protects contents and provides access to the tray handle in its vertical position.

- Made of polystyrene
- Can be secured with cable ties for extra security
- Stackable



Keeps contents clean and secure

Part No.	L x W x H	Color
HS2202	16.5 x 9.4 x 8 in (42 x 24 x 20.5 cm)	White <input type="checkbox"/>

Sample Collection Cart

Stable, secure, transport of collection supplies provides a stable, mobile center for blood collection in your facility

Standard features include:

- Made of durable stainless steel and aluminum construction
- Handle for pulling the unit
- Painted surfaces utilize Microban® protection
- 5 wheel, weighted base for stability of unit
- Designed to hold the droplet™ blood collection station

Sample Collection Cart

Part No.	HS120263
Adjustable Height	30 to 55 in (76.2 x 139.7 cm)
Designated Tray Area	16 x 9 x 6 in (42.9 x 26.5 x 17.2 cm)
Working Surface Area	14.25 x 7.5 in (36.2 x 19.1 cm)
Autoclavable Bin	—
Waste Container w/Lid	—

Deluxe Sample Collection Cart

Part No.	HS120264
Adjustable Height	30 to 55 in (76.2 x 139.7 cm)
Designated Tray Area	16 x 9 x 6 in (42.9 x 26.5 x 17.2 cm)
Working Surface Area	14.25 x 7.5 in (36.2 x 19.1 cm)
Autoclavable Bin	Clear – 15.5 x 10 x 6 in (40.6 x 25.4 x 15.2 cm) Blue – 10 x 7.4 x 6 in (25.4 x 17.8 x 15.2 cm)
Waste Container w/Lid	1.75 quarts (1.65 L)



Deluxe Sample Collection Cart

*Droplet™ blood collection tray and accessories not included