



designed for scientists



AS 4000.1 Universal attachment

/// Data Sheet

For various types of vessels. Infinitely variable clamping rolls allow universal adaptation to the vessels.

Dimensions (W x D x H): 480 x 500 x 120 mm

Set-up plate: 380 x 410 mm

Flask capacity: 54 x AS 2.1 (25 ml), 35 x AS 2.2 (50 ml), 24 x AS 2.3 (100 ml), 12 x AS 2.4 (250 ml), 9 x AS 2.5 (500 ml), 6 x AS 2.6 (1000 ml), 4 x AS 2.7 (2000 ml)

Weight: 3.200 g

Clamping screws for fixing the attachment are included in the scope of delivery of the KS 4000 incubator shaker.



designed for scientists

Technical Data

Width [mm]	55
Depth [mm]	35.1
Weight [kg]	3.375
Permissible ambient temperature [°C]	5 - 40