



designed for scientists



AS 501.2 Separating funnel attachment

/// Data Sheet

For shaking out, eluting, extracting, gassing out, dissolving, enriching, etc. Adjustment for the clamping rolls is infinitely variable, the set-up height can be changed by means of clamping devices.

Dimensions (W x D x H): 480 x 505 x 190 mm

Capacity: 12 x 50 ml or 10 x 100 ml or 6 x 250 ml.

Weight: 4180 g



designed for scientists

Technical Data

Dimensions (W x H x D) [mm]	27 x 30 x 65
Weight [kg]	4.016
Permissible ambient temperature [°C]	5 - 40