



designed for scientists



T 653 Telescopic stand

/// Data Sheet

Specially designed for the dispersing instrument T 65 digital. The stand is equipped with a pneumatic spring which enables effortless raising and lowering of the dispersion unit.

Diameter of support rod :48 mm

Dimensions (W x D): 500 x 575 mm

Height: 1.200 mm

Stroke: 500 - 1.000 mm



designed for scientists

Technical Data

Seating Plate Dimension [mm]	500 x 425
Stroke max. [mm]	1000
Dimensions (W x H x D) [mm]	500 x 1185 x 575
Weight [kg]	25.25
Permissible ambient temperature [°C]	5 - 40