



designed for scientists



R 1402 Dissolver

/// Data Sheet

The R 1402 dissolver head is designed for the use with the overhead stirrers EUROSTAR 20 high speed digital / control and the high speed stirring shaft R 50 of the disperser T 50 digital. It is ideally suited for applications which work with low to medium mixing viscosities that require radial flow patterns at high speed with medium to high shearing forces.

Working range: 1 - 30 l

Rotor diameter: 40 mm

Note: R 1402 needs to be attached to either the precision shaft R 6000 (if used with EUROSTAR 20 high speed) or high



designed for scientists

speed stirring shaft R 50 (if used with T 50 digital) in order to operate. The shafts are not included in the scope of delivery of R 1402.



designed for scientists

Technical Data

Stirrer diameter [mm]	40
Weight [kg]	0.024
Permissible ambient temperature [°C]	5 - 40