

Optional Extras Page

Double Glazed Door – available for most models

Inner Perspex Door – available for most models

Extra Shelves – available for all models

Interior illumination, Daylight or UV with 24 hour timer – all models

Interior illumination, Daylight or UV with 7 day timer – all models

Interior illumination, Daylight or UV programmable – programmable models only

Fruit Fly Modification – available for all models

Adjustable High/Low Alarm – available for all models

Access Porthole – available for all models

Integral Chart Recorder – available for all models

Digital Defrost – available for models in Series 1A & 2

Internal 13 amp socket – available for all models

Set of Castors – available for Series 4 only

RS485 Interface – available for Programmable models only

Pre-programming service – available for programmable models only

Configuration Software – available for programmable models only

Whilst the above list is not excessive we would welcome the opportunity to discuss your individual needs and assist wherever possible.