


1. Part Number: 1030A all sizes


Description

Tips for use with 12.5mm diameter PowderThief Bodies

Photo	Part No:	Nominal Volume	Weight
	1030AT-025	0.25ml	18g
	1030AT-050	0.50ml	19g
	1030AT-075	0.75ml	20g
	1030AT-100	1.00ml	24g
	1030AT-150	1.50ml	26g
	1030AT-200	2.00ml	31g

Materials of Construction:

Tip	316L stainless steel
Collar	PTFE (Part No. SS002000)


Finish: <0.5 microns Ra

BSE/TSE statement: All polishing compounds used in the manufacture of this product are free of animal products

Recommended Storage Conditions: Dry and ambient temperature

2. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Nov-20