

**1. Part Number:** 5005C-100

**2. Photo:**



### 3. Product Information:

Part No.	5005C-100
Description	Stainless Steel Dipper – 100ml
Material of Construction:	316L stainless steel
Method of Construction:	Fully Welded Construction
Surface Finish:	Cup better than 0.5 microns Ra Handle better than 1.0 microns Ra
Nominal Cup Volume:	100ml
Brimful Cup Volume:	100ml
Overall Length:	1500mm
Max. Width (cup and handle)	49mm
Cup Diameter:	39mm
Cup Height:	105mm
Nominal Weight:	1070g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-Nov-19