CASSETTES



Qty/Cs Cat. # Cat. # Color M507-2 M508-2 White M507-3 M508-3 M507-4 M508-4 1000 Green M507-5 M508-5 Yellow M507-6 M508-6 Blue M507-7 M508-7 Peach M507-8 M508-8 M507-9 M508-9 M507-10 M508-10 M507-11 M508-11 Orange M507-12 M508-12

Fluo Pink. Fluo Green and Fluo Yellow are also available. Minimum quantities apply. Contact us for more details.

M507

MICROMESHTM

Biopsy Processing / Embedding Cassettes

Made of acetal

This version of the Micromesh™ offers 1676 square openings (0.38 mm) allowing for a greatly improved fluid exchange and drainage without having to use biopsy pads. Large anterior slots in both cassette and cover ensure that the Micromesh™ biopsy cassette will sink rapidly. A large square compartment measuring 27 mm is perfect even for needle biopsies. The cover does not protrude above the cassette, a great space saving feature allowing more cassettes to be stacked in automatic labeling machines and tissue processors.

Molded from a special high density polymer, these patented cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. They can be opened and closed as often as necessary and they always relock securely without danger of specimen loss. Anterior writing area is at a 45° angle to make the cassette more suitable to be used with automated cassette printers. Available in 11 colors.

Each case contains four dispenser boxes of 250 cassettes with covers assembled.



Anatomy of Micromesh™ Cassettes

A recessed cover is a great space saving feature, allowing more cassettes to be stacked in automatic labeling machines and storage drawers



- · · · · Cover pre-installed on cassette
 - Large anterior and posterior slots in both cassette and cover ensure that the Micromesh™ biopsy cassette will be submerged rapidly

NO BIOPSY PADS REQUIRED

Anterior writing area is at a 45° angle

M508

MICROMESHTM

Biopsy Processing / Embedding Cassettes

Made of acetal

This model is similar to Series M507 but cassettes have four square compartments each measuring 13 mm. Cover and base have about 1676 square openings maximizing fluid exchange and ensuring proper drainage.

Each case contains four dispenser boxes of 250 cassettes with covers assembled.





Compatible with all cassette printers

