

SLS Lab Pro 20T3 PurA-Q3

Complete All-In-One Type 3 Water System

- Reverse Osmosis (RO)
- <75µS/cm water quality
- Up to 98% rejection of incoming feed inorganics
- Space saving footprint with integral 35L tank
- Make up rate 20L/hr
- Low energy consumption <40w
- Quiet internal boost pump and patented anti-vibration mountings
- Mercury-free option available
- First consumables included with unit
- Simple cartridge replacement



Code	Description
WAT2406	SLS Lab Pro 20T3 (PurA-Q3)
WAT2419	SLS Lab Pro 20T3 (PurA-Q3) with UV kit

Accessories

Code	Description
WAT2452	Wall mounting kit

Consumables

Code	Description
WAT2436	SLS Lab Pro pre-treatment module (RO) 20L/hr – *annual change
WAT2440	SLS Lab Pro air-vent filter (tank) – annual change
WAT2442	SLS Lab Pro UV lamp consumable for WAT2419 – annual change

*Estimate only, depends on water usage and quality of feed water.



Applications

Autoclave feed
Glasswasher feed
Feed to ultrapure
Hydroponics
Steam generators

Model Specifications

Model	DI Make-up Rate	Inorganics Rejection Rate	Organics Rejection Rate	Bacteria Rejection Rate	pH	Daily Usage (Max)
SLS Lab Pro 20T3	20L/hr	98%	>99%	99%	Neutral	200L/Day

Dimensions			Feed Water Requirements		Services Required
Width	Height	Depth	Pressure	- 5 Bar	
485mm	500mm	330mm	Conductivity	<1400 µS-cm	
Weight (Empty)			Free-Chlorine & Iron	0.1 Mg/L	
20kg			Colloid Index	<3 SDI	

00160