



Product Data Sheet

SLS8106/SLS8110

General Descriptions

Product Name	SLS Select centrifuge tube, 50mL
Further Description	Conical bottom, polypropylene (PP), white printed graduations with flat cap
General Applications	Manufactured from virgin ultra clear polypropylene to ensure good visibility Leak-resistant seal The cap grip assures easy opening and closing when wearing gloves White printed graduations with large marking strip for labelling
Container Material	Polypropylene
Cap Material	HDPE
Packaging	SLS8106 - Bagged SLS8110 - Recyclable card racks

General Features

Graduation Range	5 - 50mL
Graduation Interval	2.5mL
Temperature	-80 to 121°C
DNase/RNase	Free
Non-Pyrogenic	<0.25EU/mL
Solvent Resistant Printing	Will not rub off
Leak-Proof at	-95kPa (14psi)
Max RCF	12,000xg Refers to liquids with a specific gravity of 1.0, using appropriate rotors, cushions, and adapters

