



PRODUCT SPECIFICATION

DUAL LAYER POLYESTER ISOLATOR BONNET STERILE



PRODUCT INFORMATION

PRODUCT DESCRIPTION

Universal fit isolator cleaning tool made from dual layer continuous filament polyester wipes designed specifically to help clean inside of isolators and LAF cabinets.

PRODUCT FEATURES

Low linting, dual layer 100% continuous filament polyester for increased absorbency; Universal fit; Supplied sterile one per bag ready for use; Aqueous processing at ISO 14644-1 EU GMP Grade B; Tested to ASTM-D-629.

MANUFACTURING

Aqueous processing at ISO 14644-1 EU GMP Grade B; Packaging is performed within an ISO 14644-1 Class 7 (in use) cleanroom.

PACKAGING

Primary packaging: Double bagged in polythene packaging; Supplied 40 individually packed bonnets per carton.

STERILISATION DETAILS

Gamma Irradiation (minimum 25kGy)

QUALITY ASSURANCE

Product supplied under BS EN ISO 9001 registered systems.

SHELF LIFE

2 years from date of manufacture; Use on day of opening.

STORAGE CONDITIONS

Store in cool, well ventilated area.

MATERIAL DIMENSIONS

Wipe Material: 100% continuous filament polyester.
Wipe Weight: 250g/m²

ORDERING/PACKING INFORMATION

Sterile Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
ISED0140	Dry Dual Layer Polyester Isolator Bonnet Sterile	One Size	Each	40	245 x 195 x 260mm	1 Kg