




PRODUCT SPECIFICATION

PVC OVERSHOES NON-STERILE

TuffGuard **3**

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Disposable non-vented overshoe

PRODUCT FEATURES

Durable
Cost effective
14" or 16" sizing

PACKAGING

Primary packaged in polythene packaging

QUALITY ASSURANCE

Product supplied under BS EN ISO 9001 registered systems

STORAGE CONDITIONS

Store in a dry, well-ventilated area in ambient conditions

MATERIAL AND PERFORMANCE

Red PVC Overshoes

Polyvinylchloride (PVC) material
DIN53505 tensile strength $\geq 16.0\text{MPa}$ (MD); $\geq 15.0\text{MPa}$ (CD)
DIN53363 tear strength $\geq 85\text{N/mm}$ (MD); $\geq 117\text{N/mm}$ (CD)
DIN53455 elongation at break $\geq 190\%$ (MD); $\geq 220\%$ (CD)

White PVC Overshoes

Polyvinylchloride (PVC) material
DIN53505 tensile strength $\geq 16.8\text{MPa}$ (MD); $\geq 15.0\text{MPa}$ (CD)
DIN53363 tear strength $\geq 85\text{N/mm}$ (MD); $\geq 117\text{N/mm}$ (CD)
DIN53455 elongation at break $\geq 190\%$ (MD); $\geq 220\%$ (CD)

ORDERING/PACKING INFORMATION

Non-Sterile Code	Sterile Code	Description	Size	Pack Qty	Case Qty
SNWH81006		PVC 6"x14" Overshoe White Non-sterile	6 x 14"	500	500
SNWH83016		PVC 6"x16" Overshoe White Non-sterile	6 x 16"	500	500
SNWH82009		PVC 9"x14" Overshoe White Non-sterile	9 x 14"	500	500
SNWH81016		PVC 9"x16" Overshoe White Non-sterile	9 x 16"	500	500
SNRE81009		PVC 9"x16" Overshoe Red Non-sterile	9 x 16"	500	500